



BD™ sharps collectors

Product catalogue

Delivered **clean & new** every time



Small volume

Ideal for use in:

- Blood collection trays
- Laboratories
- Physician and dental offices
- Community



Catalogue no.	Volume (litres) and description	Dimensions H x W x D (mm)	Case quantity
302434	0.45 L	170 x 105 x 50	100
305008	0.87 L	165 x 85 x 70	50
303511	1.4 L High top	230 x 115 x 105	36
303205	1.4 L One-piece	190 x 90 x 90	36
305011	1.4 L Temporary closure	210 x 80 x 80	36
305624	1.5 L	205 x 165 x 160	40

Compatible brackets

Desktop bracket for 302434	Small cup bracket for 303511 & 305008	Large cup bracket for 303205 & 305011	Lockable cabinet for 303205 & 305011
----------------------------	---------------------------------------	---------------------------------------	--------------------------------------



305964



305970



305451



305708

Desktop bracket for 305624	Hanging bracket for 303205 & 305011	Hanging bracket for 305624
----------------------------	-------------------------------------	----------------------------



305645



305697



305710

Nestable

- Available in Fluted top or Open top options
- Fluted top is designed to limit accessibility to the collector contents
- Open top features a larger diameter to accommodate disposal needs in non-patient or non-public areas.



Catalogue no.	Volume (litres) and description	Dimensions H x W x D (mm)	Case quantity
303506	7.6 L Fluted top	250 x 280 x 165	24
303507	13.2 L Fluted top	340 x 275 x 220	20
303508	13.2 L Open top	340 x 275 x 220	20
303509	22.7 L Open top	430 x 325 x 217	12
300740	22.7 L Fluted top	430 x 325 x 217	12

Compatible brackets

Wall mount bracket for 300740, 303507, 303508 & 303509	Mobile trolley for 300740, 303507, 303508 & 303509	Hanging bracket for 303506
		
305485	305702	305706

One-piece range



- Single piece moulding, requires no assembly
- Fluted opening designed to limit accessibility to the collector contents
- Multiple bracketing options for stability

One-piece range

Catalogue no.	Volume (litres) and description	Dimensions H x W x D (mm)	Case quantity
305011	1.4 L Temporary Closure	210 x 80 x 80	36
303205	1.4 L One-piece	190 x 90 x 90	36
303207	3 L One-piece	139 x 212 x 330	24
303208	8 L One-piece	343 x 232 x 134	12
303209	17 L One-piece	450 x 266 x 190	8

Compatible brackets

Cup bracket for 303205 & 305011	Quick release strap for 303207, 303208 & 303209	Hanging bracket for 305011 & 303205
---------------------------------	---	-------------------------------------



305451



305497



305697

Hanging bracket for 303207 & 303208	Wall mount bracket for 303207 & 303208	Hanging bracket for 303209
-------------------------------------	--	----------------------------



305698



305699



305700

Wall mount bracket for 303209	Mobile trolley for 303209	Lockable cabinet for 305011 & 303205
-------------------------------	---------------------------	--------------------------------------



305701



305702



305708

Counter Balance & Next Generation

- Designed to improve tamper resistance by preventing hand access to the contents
- Compatible with locking cabinets to provide security and prevent unwanted access to collector contents



Catalogue no.	Volume (litres) and description	Dimensions H x W x D (mm)	Case quantity
301272	5.1 L Next Generation	272 x 297 x 108	20
301273	7.6 L Next Generation	289 x 333 x 150	10
301274	11.3 L Next Generation	420 x 333 x 150	10



Catalogue no.	Volume (litres) and description	Dimensions H x W x D (mm)	Case quantity
305081	7.6 L Counter Balance	289 x 333 x 150	10

Counter Balance & Next Generation

Compatible brackets

Small locking cabinet for 301272	Small BD Recykleen™ locking cabinet for 301272	BD Recykleen™ locking cabinet for 301273 & 301274
----------------------------------	--	---



305017



305096



305097

Wall mount bracket for 301273, 305081 & 301274	Locking cabinet for 301273 & 301274	Wall mount bracket for 301272
--	-------------------------------------	-------------------------------



305409



305420



305447

Mobile trolley for 301273, 305081 & 301274	Hanging bracket for 301272
--	----------------------------



305702



305703

Mobile point of use unit for 301272	Hanging bracket for 301273, 305081 & 301274
-------------------------------------	---



305704



305705

BD Recykleen™ cabinets are made from 100% recycled polycarbonate plastic.

Large capacity

Ideal for use in:

- Operating rooms
- Laboratories
- Dialysis units



Catalogue no.	Volume (litres) and description	Dimensions H x W x D (mm)	Case quantity
305608	34 L Slide top	479 x 447 x 298	8
302406	72 L Slide top	660 x 488 x 357	5

Compatible BD Recykleen™ trolleys

BD Recykleen™ foot-operated trolley 34 L for 305608

BD Recykleen™ foot-operated trolley 72 L for 302406



305091



305093

BD Recykleen™ trolleys

- Designed to help minimise the risk of collectors tipping over
- Wheels enable easy portability of the collector
- Locking wheels, ensuring the product is mobile only when desired
- Designed with the ability to open and close the collector for hands-free disposal

Chemotherapy



Catalogue no.	Volume (litres) and description	Dimensions H x W x D (mm)	Case quantity
303210	8 L One-piece	343 x 232 x 134	12
303510	17 L One-piece	450 x 266 x 190	8

Compatible brackets

Hanging bracket for 303210	Wall mount bracket for 303210	Quick release strap for 303210 & 303510
----------------------------	-------------------------------	---



305698



305699



305497

Wall mount bracket for 303510	Hanging bracket for 303510	Mobile trolley for 303510
-------------------------------	----------------------------	---------------------------



305700



305701



305702

BD sharps collector dimensions

Style	Size	Catalogue No.	Sharps Collector Dimensions, H x W x D (mm)
Small volume	0.45 L	302434	170 x 105 x 50
	0.87 L	305008	165 x 85 x 70
	1.4 L High top	303511	230 x 115 x 105
	1.4 L	303205	190 x 90 x 90
	1.4 L Temp closure	305011	210 x 80 x 80
	1.5 L	305624	205 x 165 x 160
Nestable	7.6 L Fluted top	303506	250 x 280 x 165
	13.2 L Fluted top	303507	340 x 275 x 220
	13.2 L Open top	303508	340 x 275 x 220
	22.7 L Fluted top	300740	430 x 325 x 217
	22.7 L Open top	303509	430 x 325 x 217
One-piece	1.4 L	303205	190 x 90 x 90
	1.4 L Temp closure	305011	210 x 80 x 80
	3 L	303207	139 x 212 x 330
	8 L	303208	343 x 232 x 134
	17 L	303209	450 x 266 x 190
Counter Balance	7.6 L	305081	289 x 333 x 150
Next Generation	5.1 L	301272	272 x 297 x 108
	7.6 L	301273	289 x 333 x 150
	11.3 L	301274	420 x 333 x 150
Large capacity	34 L Slide top	305608	479 x 447 x 298
	72 L Slide top	302406	660 x 488 x 357
Chemotherapy	8 L One-piece	303210	343 x 232 x 134
	17 L One-piece	303510	450 x 266 x 190

BD brackets with compatible collectors

Catalogue No.	Description	Sale UOM	Compatible collector
305017	Small locking cabinet	EA	301272
305091	BD Recykleen™ foot-operated trolley 34 L	EA	305608
305093	BD Recykleen™ foot-operated trolley 72 L	EA	302406
305096	BD Recykleen™ small locking cabinet	EA	301272
305097	BD Recykleen™ locking cabinet	EA	301273 & 301274
305409	Wall mount bracket	EA	301273, 301274 & 305081
305420	Locking cabinet	EA	301273 & 301274
305447	Wall mount bracket	EA	301272
305451	Large c up bracket	EA	303205 & 305011
305485	Wall mount bracket	EA	300740 303507, 303508 & 303509
305497	Quick release strap	EA	303207, 303208, 303209, 303210 & 303510
305645	Desktop bracket	EA	305624
305697	Hanging bracket	EA	303205 & 305011
305698	Hanging bracket	EA	303207, 303208 & 303210
305699	Wall mount bracket	EA	303207, 303208 & 303210
305700	Hanging bracket	EA	303209
305701	Wall mount bracket	EA	303209 & 303510
305702	Mobile trolley	EA	303507, 303508, 303509, 303510 300740, 301273, 305081 & 301274
305703	Hanging bracket	EA	301272
305704	Mobile point of use unit	EA	301272
305705	Hanging bracket	EA	301273, 301274 & 305081
305706	Hanging bracket	EA	303506
305708	Lockable cabinet	EA	303205 & 305011
305710	Hanging bracket	EA	305624
305964	Desktop bracket	EA	302434
305970	Small cup bracket	EA	303511 & 305008

Safe sharps disposal tips

- Do not fill past the fill line
- Change the sharps collector as soon as the contents reach the fill line
- When removing a sharps collector, always use the handles provided
- When assembly is required, secure the lid to the sharps collector base correctly, an audible click will be heard once the lid is in place
- Do not dispose of non-sharps waste in sharps collectors as this adds unnecessary cost to your waste disposal fees
- Always remember to aim for point of use positioning
- When placing a sharp into a sharps collector, place the sharp end in first. If using a Counter Balance or Next Generation sharps collector, you will still place the sharp end in first, but on a horizontal angle
- Do not recap needles
- Do not de-notch needles off syringes
- If you see that a collector has been positioned unsafely please alert the relevant person in your department